

Recombinant human 39S ribosomal protein L11, mitochondrial

Catalog No: #AP71632

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP71632-1 20ug #AP71632-2 100ug #AP71632-3 1mg

Description

Product Name	Recombinant human 39S ribosomal protein L11, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-191aaSequence Info:Partial
Accession No.	Q9Y3B7
Uniprot	Q9Y3B7
GeneID	65003;
Calculated MW	36.6 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MSKLGRAARGLRKPEVGGVIRAIVRAGLAMPGPPLGPNVLRGQVRSINQFCKEFNERTKDIKEGIPLPTKILVKP DRTFEIKIGQPTVSYFLKAAAGIEKGARQTGKEVAGLVTLKHVYEIARIKAQDEAFALQDVLSSVVRSSIQSARS LGIRVVKDLSSEELAAFQKERAIFLAAQKEADLAAQEAAK
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Structural compensation for the deficit of rRNA with proteins in the mammalian mitochondrial ribosome. Systematic analysis of protein components of the large ribosomal subunit from mammalian mitochondria.Suzuki T., Terasaki M., Takemoto-Hori C., Hanada T., Ueda T., Wada A., Watanabe K.J. Biol. Chem. 276:21724-21736(2001)Research Topic:Epigenetics and Nuclear Signaling

Note: This product is for in vitro research use only