

Recombinant E.coli 30S ribosomal protein S2

Catalog No: #AP72026



Package Size: #AP72026-1 20ug #AP72026-2 100ug #AP72026-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E.coli 30S ribosomal protein S2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:2-241aaSequence Info:Full Length
Accession No.	P0A7V0
Uniprot	P0A7V0
GeneID	947874;
Calculated MW	53.6 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	ATVSMRDMLKAGVHFGHQTRYWNPKMKPFIFGARNKVHIINLEKTVPMFNEALAE LNKI ASRK GKILFVGTKRA ASEAVKDAALSCDQFFVNHRWLGGMLTNWKTVRQSIKRLKDLQSQDGTDFDKLTKKEALMRTRELEKLENS LGGIKDMGGLPDALFVIDADHEHIAIKEANNL GIPVFAIVDTNSDPDGVD FVIPGNDDAIRAVTLYLGAVAATVRE GRSQDLASQAEE SFVEAE
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Highly accurate genome sequences of Escherichia coli K-12 strains MG1655 and W3110.Hayashi K., Morooka N., Yamamoto Y., Fujita K., Isono K., Choi S., Ohtsubo E., Baba T., Wanner B.L., Mori H., Horiuchi T.Mol. Syst. Biol. 2:E1-E5(2006)Research Topic:Others

Note: This product is for in vitro research use only