

Recombinant E.coli 50S ribosomal protein L33

Catalog No: #AP72040



Package Size: #AP72040-1 20ug #AP72040-2 100ug #AP72040-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E.coli 50S ribosomal protein L33
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:2-54aaSequence Info:Partial
Accession No.	P0A7N9
Uniprot	P0A7N9
GeneID	8913847;946318;
Calculated MW	10.1 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	AKGIREKIKLVSSAGTGHFYTTTKNKRTKPEKLELKKFDPVVRQHVIIYKEAKI
Formulation	Tris-based buffer50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.</p>

References

Observation of Escherichia coli ribosomal proteins and their posttranslational modifications by mass spectrometry.Arnold R.J., Reilly J.P.Anal. Biochem. 269:105-112(1999)Research Topic:Others

Note: This product is for in vitro research use only