

Recombinant E.coli 50S ribosomal protein L28

Catalog No: #AP72047



Package Size: #AP72047-1 20ug #AP72047-2 100ug #AP72047-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E.coli 50S ribosomal protein L28
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:2-78aaSequence Info:Full Length
Accession No.	P0A7M2
Uniprot	P0A7M2
GeneID	946941;
Calculated MW	35.9 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	SRVCQVTGKRPVTGNNRSHALNATKRRFLPNLHSHRFWVESEKRFVTLRVSAKGMRVIDKKGIDTVLAELRA RGEKY
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Observation of Escherichia coli ribosomal proteins and their posttranslational modifications by mass spectrometry.Arnold R.J., Reilly J.P.Anal. Biochem. 269:105-112(1999)Research Topic:Others

Note: This product is for in vitro research use only