

Recombinant Gallus gallus Uncharacterized protein

Catalog No: #AP72282



Package Size: #AP72282-1 20ug #AP72282-2 100ug #AP72282-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Gallus gallus Uncharacterized protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-361aaSequence Info:Full Length
Accession No.	Q98SS5
Uniprot	Q98SS5
Calculated MW	44 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MELDFGHFDERDKASRNMRGTRMNGLPSPTHSAHCSFYRTRTLQALSNEKKAKKVRFYRNGDRYFKGIVYA VSSDRFRSFDALLADLTRSLSDNINLPQGVRYIYIDGSRKIGSMDELEEGESYVCSSDNFFKKVEYTKNVNPN WSVNVKTSANQKAPQSLASSNSAQAKENKDFVRPKLVTIIRSGVKPRKAVRVLLNKKTAHSFEQVLTDITEAIK LETGVVKKLYTLDGKQVTCLHDFFGDDVFIACGPEKFRYAQDDFSLDESECRMKGSPAAATGSKSSPTPQ KSSAKSPAPMRRSKSPADSANGTSSSQLSTPKSKQSPISTPTSPGSLRKHKVDLYLPLSLDDSDSLGDSM
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Screening from a subtracted embryonic chick hindbrain cDNA library identification of genes expressed during hindbrain, midbrain and cranial neural crest development.Christiansen J.H., Coles E.G., Robinson V., Pasini A., Wilkinson D.G.Mech. Dev. 102:119-133(2001)Research Topic:Others

Note: This product is for in vitro research use only