

Recombinant rat Interferon gamma protein

Catalog No: #AP72295



Package Size: #AP72295-1 20ug #AP72295-2 100ug #AP72295-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant rat Interferon gamma protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:23-156aaSequence Info:Full Length
Accession No.	P01581
Uniprot	P01581
GeneID	25712;
Calculated MW	42.5 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	QGTLIESLES LKNYFNSSMDAMEGKSLLLDIWRNWQKDGNTKILESQIISFYLR LFEVLKDNQAISNNISVIESH LITNFFSNSKAKKDAFMSIAKFEVNNPQIQHKAVNELIRVIHQLSPESSLRKRKR SRC
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Produced by lymphocytes activated by specific antigens or mitogens. IFN-gamma, in addition to having antiviral activity, has important immunoregulatory functions. It is a potent activator of macrophages, it has antiproliferative effects on transformed cells and it can potentiate the antiviral and antitumor effects of the type I interferons.

References

Cloning and expression of the chromosomal immune interferon gene of the rat.Dijkema R., van der Meide P.H., Pouwels P.H., Caspers M., Dubbeld M., Schellekens H.EMBO J. 4:761-767(1985)Research Topic:Others

Note: This product is for in vitro research use only