

## Recombinant human BoIA-like protein 1

Catalog No: #AP72458



Package Size: #AP72458-1 20ug #AP72458-2 100ug #AP72458-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant human BoIA-like protein 1
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-137aaSequence Info:Full Length
Other Names	hBoIA
Accession No.	Q9Y3E2
Uniprot	Q9Y3E2
GeneID	51027;
Calculated MW	16.3 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MLSGRLVLGLVSMAGRVCCLCQGSAGSGAIGPVEAAIRTKLEEALSPEVLELRNESGGHAVPPGSETHFRVAVV SSRFEGLSPLQRHRLVHAALAEELGGPVHALAIQARTPAQWRENSQLDTSPPCLGGNKKTGLTP
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## References

Identification of novel human genes evolutionarily conserved in *Caenorhabditis elegans* by comparative proteomics.Lai C.-H., Chou C.-Y., Ch'ang L.-Y., Liu C.-S., Lin W.-C.Genome Res. 10:703-713(2000)Research Topic:Cell Biology

Note: This product is for in vitro research use only