

Recombinant Homo sapiens Lutropin-choriogonadotropic hormone receptor

Catalog No: #AP72649

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP72649-1 20ug #AP72649-2 100ug #AP72649-3 1mg

Description

Product Name	Recombinant Homo sapiens Lutropin-choriogonadotropic hormone receptor
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:27-363aaSequence Info:Partial
Other Names	Luteinizing hormone receptor ;LHR ;LSH-R
Accession No.	P22888
Uniprot	P22888
GeneID	3973;
Calculated MW	39.6 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	EALCPEPCNCVPD GALRCPGPTAGLTRLRLSLAYLPVKVIPSQAFRGLNEVIKIEISQIDSLERIEANAFDNLNLSEI LIQNTKNLRYIEPGA FINLPRLKYLSICNTGIRKFPDVTKVFSSES NFILEICDNLHITIPGNAFQGMNNE SVTLKL YGNGFEEVQSHAFNGTTLT SLELKENVHLEKMHNGAFRGATGPKTLDISSTKLQALPSYGLESIQR LIATSSYS LKKLPSRETFVNLEATLTYP SHCCAFRNLP TKEQNF SHSISENF SKQCESTVRKVNNK TLYSSMLAESEL SG WDYEYGFCLPKTPRCAPEPDAFNPCEDIMGYDFLR
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Receptor for lutropin-choriogonadotropic hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

References

Cloning and sequencing of human LH,hCG receptor cDNA.Minegish T., Nakamura K., Takakura Y., Miyamoto K., Hasegawa Y., Ibuki Y., Igarashi M.Biochem. Biophys. Res. Commun. 172:1049-1054(1990)Research Topic:Neuroscience

Note: This product is for in vitro research use only