

Recombinant mouse Mesothelin

Catalog No: #AP72686



Package Size: #AP72686-1 20ug #AP72686-2 100ug #AP72686-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant mouse Mesothelin
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:298-600aaSequence Info:Partial
Other Names	Pre-pro-megakaryocyte-potentiating factor
Accession No.	Q61468
Uniprot	Q61468
GeneID	56047;
Calculated MW	36.1 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	DAEQKACPPGKEPYKVEDLIFYQNWLEACVDGTMRLARQMDLVNEIPFTYEQLSIFKHKLDKTYPPQGYPELSL QQLGHFFRYVSPEDIHQWNVTSPTDVKTLKLVSKGQKMNAQAIALVACYLRGGGQLDEDMVKALGDIPLSYL CDFSPQDLHSPSSVMWLVGPQDLKCSQRHLGLLYQKACSAFQNVSGLEYFEKIKTFLGGASVKDLRALSQ HNVSMDIATFKRLQVDSLVLGSLVAEVQKLLGPNIVDLKTEEDKSPVRDWLFRQHQQKDLDRGLGLQGGIPNGY LVLDFNVREAFS
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Mbrane-anchored forms may play a role in cellular adhesion.

References

Stromelysin-1 and mesothelin are differentially regulated by Wnt-5a and Wnt-1 in C57mg mouse mammary epithelial cells.Prieve M.G., Moon R.T.BMC Dev. Biol. 3:2-2(2003)Research Topic:Others

Note: This product is for in vitro research use only