

Recombinant mouse Prostate stem cell antigen

Catalog No: #AP72740



Package Size: #AP72740-1 20ug #AP72740-2 100ug #AP72740-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant mouse Prostate stem cell antigen
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:21-95aaSequence Info:Full Length
Accession No.	P57096
Uniprot	P57096
GeneID	72373;
Calculated MW	10.4 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	LQCYSTAQMNNRDCLNVQNCSLDQHSCFTSRIRAIGLVTVISKGCSSQCEDDSENYLGGKKNITCCYSDLCN VN
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May be involved in the regulation of cell proliferation.

References

Prostate stem cell antigen a cell surface marker overexpressed in prostate cancer.Reiter R.E., Gu Z., Watabe T., Thomas G., Szigeti K., Davis E., Wahl M., Nisitani S., Yamashiro J., le Beau M.M., Losa M., Witte O.N.Proc. Natl. Acad. Sci. U.S.A. 95:1735-1740(1998)Research Topic:Others

Note: This product is for in vitro research use only