

# Recombinant Rabbit Vascular endothelial growth factor A



Catalog No: #AP72840

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP72840-1 20ug #AP72840-2 100ug #AP72840-3 1mg

## Description

Product Name	Recombinant Rabbit Vascular endothelial growth factor A
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:51-511aaSequence Info:Full Length
Accession No.	XP_002714743.1
Calculated MW	51.7 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	YLWLLHAPLLHVSLDSLVAAFSVVFEVLPSSLDIPGSKKEEQASTQLGGVEEEHEASGVVTRQLAFVVDITNP TLEKETKTITRIHKAPKFP SHFGFWKRAEEERGGKSSGEKSRKTDGVGERAQAGEQREGPGQRAAALDRQT DTAPSPSAHLLPGRRPTVDAASGGQQPEPAPEGGVEGVGARGIALKLFVQLLGCSRSGGVVVRAGGAEPS GAARSVSSGREPPPPPEEEEGEEEGEKEEERGPRWRLGAGEPGSWTGAAVVCADSAPAARAPQALARAS APGARGARGGAEESGPRSPRRGSASRAGPGRASET MNFLLSWVHWSLALLLYLHHAKWSQAAPMAEEG DNKPHEVVKFMEVYRRSYCQPIETLVDIFQEY PDEIEYIFKPCVPLVRCGGCCNDESLECVPT EEFNVTMQIM RIKPHQGQHIGEMSFLQHNKCECRCDKPRR
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## References

Research Topic:Others

Note: This product is for in vitro research use only