

Recombinant Putative defensin-like protein 70

Catalog No: #AP72861



Package Size: #AP72861-1 20ug #AP72861-2 100ug #AP72861-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Putative defensin-like protein 70
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:28-82aaSequence Info:Full Length
Other Names	Putative low-molecular-weight cysteine-rich protein 83 Short name: Protein LCR83
Accession No.	P82792
Uniprot	P82792
GeneID	3771507;
Calculated MW	32.9 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	NFASGEASSQLCFNPCTPQLGNNECNTICMNKKYKEGSCVGFIPPTSKYCCCKT
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

"Genome organization of more than 300 defensin-like genes in Arabidopsis."Silverstein K.A.T., Graham M.A., Paape T.D., VandenBosch K.A. Plant Physiol. 138:600-610(2005)Research Topic:Others

Note: This product is for in vitro research use only