

Recombinant Homo sapiens Fibroleukin

Catalog No: #AP73048

Package Size: #AP73048-1 20ug #AP73048-2 100ug #AP73048-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Homo sapiens Fibroleukin
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:24-439aaSequence Info:Full Length
Other Names	Fibrinogen-like protein 2 pT49
Accession No.	Q14314
Uniprot	Q14314
GeneID	10875;
Calculated MW	49.6 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	<p>NNETEEIKDERAKDVCVRLERLSESRGKCEEAGECPYQVSLPPLTIQLPKQFSRIEEVFKEVQNLKEIVNSLKKSCQ DCKLQADDNGDPGRNGLLLPSTGAPGEVGDNRVRELESEVNLKSELKNAKEEINVLHGRLEKLNLVNMNIE NYVDSKVANLTFVNSLDGKCSKCPSEQIQSRPVQHLYKDCSDYYAIGKRSSSETYRVTPDPKNSSFEVYCD METMGGGWTVLQARLDGSTNFRTRWQDYKAGFGNLRREFWLGNDKIHLTKSKEMILRIDLEDFNGVELYAL YDQFYVANEFLKYRLHVGNYNGTAGDALRFNKHYNHDLKFFFTDPKDNDRYPSGNCGLYSSGWWFACL ANLNGKYYHQYRGVNRNGIFWGTWPGVSEAHPPGGYKSSFKEAKMMIRPKHFKP</p>
Formulation	Tris-based buffer50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.</p>

Background

May play a role in physiologic lymphocyte functions at mucosal sites

References

"Glycoproteomics analysis of human liver tissue by combination of multiple enzyme digestion and hydrazide chemistry."Chen R., Jiang X., Sun D., Han G., Wang F., Ye M., Wang L., Zou H.J. Proteome Res. 8:651-661(2009)Research Topic:Immunology

Note: This product is for in vitro research use only