

## Recombinant human Protein ERGIC-53-like

Catalog No: #AP73120



Package Size: #AP73120-1 20ug #AP73120-2 100ug #AP73120-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Recombinant human Protein ERGIC-53-like
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:26-462aaSequence Info:Full Length
Other Names	ERGIC53-like protein;Lectin mannose-binding 1-like ;LMAN1-like protein
Accession No.	Q9HAT1
Uniprot	Q9HAT1
GenID	79748;
Calculated MW	49.5 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	GCPPLRRFEYKLSFKGPRLALPGAGIPFWSHHGDAILGLEEVRLTPSMRNRSGAVWSRASVPFSAWEVEVQ MRVTGLGRRGAQGMVWYTRGRGHVGSVLGGLASWDGIGIFFDSPAEDTQDSPAIRVLASDGHIPSEQPGD GASQGLGSCHWDFRNRPHFRARITYWGQRLRMSLNSGLTPSDPGFCVDVGPLLLVPGGFFGVSAATGTL ADDHDVLSFLTSLSEPSPEVPPQPFLEMQLRLARQLEGLWARLGLGTREDVTPKSDSEAQGEGERLFDLE ETLGRHRRILQALRGLSKQLAQERQWKKQLGPPGQARPDGGWALDASCQIPSTPGRGGHLSMSLNKDSAK VGALLHGQWTLQLALQEMRDAAVRMAAEAQVSYLPVGIEHHFLELDHILGLLQEELRGPAAAAKAPRPPGQ PPRASSCLQ
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## References

ERGL, a novel gene related to ERGIC-53 that is highly expressed in normal and neoplastic prostate and several other tissues.Yerushalmi N., Keppler-Hafkemeyer A., Vasmatis G., Liu X.-F., Olsson P., Bera T.K., Duray P., Lee B., Pastan I.Gene 265:55-60(2001)Research Topic:Signal Transduction

Note: This product is for in vitro research use only