

Recombinant Homo sapiens Suppressor of SWI4 1 homolog

Catalog No: #AP73125

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP73125-1 20ug #AP73125-2 100ug #AP73125-3 1mg

Description

Product Name	Recombinant Homo sapiens Suppressor of SWI4 1 homolog
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-473aaSequence Info:Full Length
Other Names	Brix domain-containing protein 3;Peter Pan homolog
Accession No.	Q9NQ55
Uniprot	Q9NQ55
GeneID	56342;692312;
Calculated MW	55.2 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MGQSGRSRHRQKRARAQAQLRNLEAYAANPHSFVFTRGCTGRNIRQLSLDVRVMEPLTASRLQVRKKNLSK DCVAVAGPLGVTHFLILSKTETNVYFKLMRLPGGPTLTFQVKKYSLVRDVSLSLRRHRMHEQQFAHPPLLVLN SFGPHGMHVKLMATMFQNLFPSINVHKVNLNTIKRCLLIDYNPDSQELDFRHSYKVVVPGASRGMKLLQEKF PNMSRLQDISELLATGAGLSESEAEPDGDHNTLQPQAVAGRGNMRAQQSAVRLTEIGPRMTLQLIKVQEGVG EGKVMFHSFVSKTEELQAILEAKEKKLRLKAQRQAQAQNVQRKQEQREHRKKSLEGMKKARVGGSDDEE ASGIPSRTASLELGEDDDEQEDDDIEYFCQAVGEAPSEDLFPEAKQKRLAKSPGRKRKRWEMDRGRGLCD QKFKPKTKDKSQGAQARRGPRGASRDGGRGRGRGRPGKRVA
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May have a role in cell growth.

References

Cloning, genomic organization, and tissue distribution of human Ssf-1.Suarez-Huerta N., Boeynaems J.-M., Communi D.Biochem. Biophys. Res. Commun. 275:37-42(2000)Research Topic:Epigenetics and Nuclear Signaling

Note: This product is for in vitro research use only