

## Recombinant Human Filaggrin

Catalog No: #AP73157



Package Size: #AP73157-1 20ug #AP73157-2 100ug #AP73157-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Recombinant Human Filaggrin
Brief Description	Recombinant Protein
Host Species	Baculovirus
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:3838-4061aaSequence Info:Partial
Accession No.	P20930
Uniprot	P20930
GeneID	2312;
Calculated MW	26.8 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	DSSRHSQSGQGESAGSRRSRRQGSSVSQDSDSEAYPEDSERRSEASARNHHGSSREQSRDGSRHGSSSH RDTASHVQSSPVQSDSSTAKEHGHFSSLSQDSAYHSGIQSRGSPHSSSSYHYQSEGTERQKQGSLVWRH GSYGSADYDYGESGFRHSQHGSVSYNSNPVVKERSDICKASAFGKDHPRYATYINKDPGLCGHSSDISKQ LGFSQSQRYYYYE
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Aggregates keratin intermediate filaments and promotes disulfide-bond formation among the intermediate filaments during terminal differentiation of mammalian epidermis.

## References

The DNA sequence and biological annotation of human chromosome 1.Gregory S.G., Barlow K.F., McLay K.E., Kaul R., Swarbreck D., Dunham A., Scott C.E., Howe K.L., Woodfine K., Spencer C.C.A., Jones M.C., Gillson C., Searle S., Zhou Y., Kokocinski F., McDonald L., Evans R., Phillips K. Bentley D.R.Nature 441:315-321(2006) Research Topic:Signal Transduction

Note: This product is for in vitro research use only