

MCTS1 Conjugated Antibody

Catalog No: #C29287



Package Size: #C29287-AF350 100ul #C29287-AF405 100ul #C29287-AF488 100ul
 #C29287-AF555 100ul #C29287-AF594 100ul #C29287-AF680 100ul
 #C29287-AF750 100ul #C29287-Biotin 100ul #C29287-F647 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	MCTS1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MCTS1 (NP_054779.1).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MCTS1;MCT-1;MCT1
Accession No.	Swiss Prot:Q9ULC4GeneID:28985
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	19kDa/20kDa
SDS-PAGE MW	20kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

-

Note: This product is for in vitro research use only