

# MonkeyB PlateletB FactorB 4B (PF4)B ELISAB Kit

Catalog No: #EK20163



Package Size: #EK20163-1 48T #EK20163-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MonkeyB PlateletB FactorB 4B (PF4)B ELISAB Kit
Applications	ELISA
Species Reactivity	Monkey

## Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.122ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only