

# Human Eukaryotic peptide chain release factor GTP-binding subunit ERF3A (GSPT1) ELISA Kit



Catalog No: #EK20165

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20165-1 48T #EK20165-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Eukaryotic peptide chain release factor GTP-binding subunit ERF3A (GSPT1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:125-8000pg/ml

Sensitivity:43pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only