

Rat Thyroid Stimulating Hormone (TSH) ELISA Kit

Catalog No: #EK20193



Package Size: #EK20193-1 48T #EK20193-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Thyroid Stimulating Hormone (TSH) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.234-15ng/ml

Sensitivity:0.072ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only