

Rat Thyroid Stimulating Hormone Receptor Antidoby (TRAb) ELISA Kit



Catalog No: #EK20213

Orders: order@signalwayantibody.com

Package Size: #EK20213-1 48T #EK20213-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat Thyroid Stimulating Hormone Receptor Antidoby (TRAb) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:125-4000pg/ml
Sensitivity:43pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only