

# Rat Anti-Thyroid-Peroxidase Antibody (TPO-Ab) ELISA Kit

Catalog No: #EK20236

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20236-1 48T #EK20236-2 96T

## Description

Product Name	Rat Anti-Thyroid-Peroxidase Antibody (TPO-Ab) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:0.625-20ng/ml

Sensitivity:0.253ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only