

ITLN1 Rabbit Polyclonal Conjugated Antibody

Catalog No: #C29672



Package Size: #C29672-AF350 100ul #C29672-AF405 100ul #C29672-AF488 100ul
 #C29672-AF555 100ul #C29672-AF594 100ul #C29672-AF647 100ul
 #C29672-AF680 100ul #C29672-AF750 100ul #C29672-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	ITLN1 Rabbit Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ITLN1 (NP_060095.2).
Other Names	ITLN1;HL-1;HL1;INTL;ITLN;LFR;hIntL;omentin
Accession No.	Swiss Prot:Q8WWA0GeneID:55600
Calculated MW	34kDa
SDS-PAGE MW	37kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only