

Human Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) ELISA Kit

Catalog No: #EK20170

Orders: order@signalwayantibody.com

Package Size: #EK20170-1 48T #EK20170-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Hu

Application Details

Detect Range:62.5-4000pg/ml

Sensitivity:21.9pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only