

## Recombinant Rabbit Interleukin-1 beta(IL1B)

Catalog No: #AP89630



Package Size: #AP89630-1 20ug #AP89630-2 100ug #AP89630-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rabbit Interleukin-1 beta(IL1B)
Brief Description	Recombinant Rabbit Interleukin-1 beta(IL1B)
Target Name	His-tagged
Other Names	IL1BInterleukin-1 beta; IL-1 beta; Lymphocyte proliferation-potentiating factor
Tag Info	N-terminal His-tagged
Target Sequence	117-268 AVRS LHCRLQDAQQ KSLVLSGTYE LKALHLNAEN LNQQVVFMS FVQGEESNDK IPVALGLRGK NLYLSCVMKD DKPTLQLESV DPNRYPKKKM EKRFVFNKIE IKDKLEFESA QFPNWIYSTS QTEYMPVFLG NNSGGQLID FSMEFVSS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

&gt;85% (SDS-PAGE)

Note: This product is for in vitro research use only