

## ACTH (7-23) Polyclonal Antibody Biotin Conjugated

Catalog No: #C90001Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ACTH (7-23) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	7-23/39
Conjugates	Bio
Target Name	ACTH (7-23)
Other Names	ACTH (7-23)
Accession No.	Uniprot:P01189Gene ID:5443
Uniprot	P01189
GeneID	5443
Calculated MW	4.3
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only