

Tyrosine Hydroxylase Polyclonal Antibody AP Conjugated

Catalog No: #C90006AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine Hydroxylase Polyclonal Antibody AP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	101-165/528
Conjugates	AP
Target Name	Tyrosine Hydroxylase
Other Names	Tyrosine Hydroxylase
Accession No.	Uniprot:P07101Gene ID:7054
Uniprot	P07101
GeneID	7054
Calculated MW	60
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200

Note: This product is for in vitro research use only