

ERK1/ERK2 Polyclonal Antibody PE Conjugated

Catalog No: #C90009PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ERK1/ERK2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-358/358
Conjugates	PE
Target Name	ERK1 / ERK2
Other Names	ERK1 / ERK2
Accession No.	Uniprot:P27361Gene ID:5595
Uniprot	P27361
GeneID	5595
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only