

# GSK-3 beta mouse Monoclonal Antibody Biotin Conjugated



Catalog No: #C90010Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GSK-3 beta mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-50/420
Conjugates	Bio
Target Name	GSK-3 beta
Other Names	GSK-3 beta
Accession No.	Uniprot:P49841Gene ID:2932
Uniprot	P49841
GeneID	2932
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only