

## Osteopontin Polyclonal Antibody Biotin Conjugated

Catalog No: #C90013Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Osteopontin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Ms Rat
Immunogen Type	221-280/294
Conjugates	Bio
Target Name	Osteopontin
Other Names	Osteopontin
Accession No.	Uniprot:P10923Gene ID:20750
Uniprot	P10923
GeneID	20750
Calculated MW	34
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only