

Melatonin Receptor 1A Polyclonal Antibody PE Conjugated

Catalog No: #C90014PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melatonin Receptor 1A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-280/350
Conjugates	PE
Target Name	Melatonin Receptor 1A
Other Names	Melatonin Receptor 1A
Accession No.	Uniprot:P48039Gene ID:4543
Uniprot	P48039
GeneID	4543
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only