

Bcl-2 Polyclonal Antibody HRP Conjugated

Catalog No: #C90016HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bcl-2 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	101-160/236
Conjugates	HRP
Target Name	Bcl-2
Other Names	Bcl-2
Accession No.	Uniprot:P49950Gene ID:596
Uniprot	P49950
GeneID	596
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200

Note: This product is for in vitro research use only