

## AIF mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C90020FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AIF mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	131-230/613
Conjugates	FITC
Target Name	AIF
Other Names	AIF
Accession No.	Uniprot:O95831Gene ID:9131
Uniprot	O95831
GeneID	9131
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only