

Insulin Receptor alpha Polyclonal Antibody PE Conjugated

Catalog No: #C90023PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin Receptor alpha Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	701-760/1382
Conjugates	PE
Target Name	Insulin Receptor alpha
Other Names	Insulin Receptor alpha
Accession No.	Uniprot:P06213Gene ID:3643
Uniprot	P06213
GeneID	3643
Calculated MW	80/152
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only