

## Melan A Polyclonal Antibody Biotin Conjugated

Catalog No: #C90026Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Melan A Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu
Immunogen Type	51-118/118
Conjugates	Bio
Target Name	Melan A
Other Names	Melan A
Accession No.	Uniprot:Q16655Gene ID:2315
Uniprot	Q16655
GeneID	2315
Calculated MW	13
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only