

CALCA Polyclonal Antibody PE Conjugated

Catalog No: #C90033PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CALCA Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Rat
Immunogen Type	121..141/141
Conjugates	PE
Target Name	CALCA
Other Names	CALCA
Accession No.	Uniprot:P01258Gene ID:796
Uniprot	P01258
GeneID	796
Calculated MW	44240
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only