

Neurokinin 1 Receptor Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90034PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Neurokinin 1 Receptor Polyclonal Antibody PE-Cy7 Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 2-50/407 |
| Conjugates | PE-Cy7 |
| Target Name | Neurokinin 1 Receptor |
| Other Names | Neurokinin 1 Receptor |
| Accession No. | Uniprot:P25103Gene ID:6869 |
| Uniprot | P25103 |
| GeneID | 6869 |
| Calculated MW | 45 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Note: This product is for in vitro research use only