

Substance P Polyclonal Antibody Biotin Conjugated

Catalog No: #C90035Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Substance P Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	58-68/129
Conjugates	Bio
Target Name	Substance P
Other Names	Substance P
Accession No.	Uniprot:P20366Gene ID:6863
Uniprot	P20366
GeneID	6863
Calculated MW	1.4/13
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only