

TNF alpha Polyclonal Antibody HRP Conjugated

Catalog No: #C90042HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | TNF alpha Polyclonal Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IHC |
| Species Reactivity | Hu |
| Immunogen Type | 86-150/233 |
| Conjugates | HRP |
| Target Name | TNF alpha |
| Other Names | TNF alpha |
| Accession No. | Uniprot:P01375Gene ID:7124 |
| Uniprot | P01375 |
| GeneID | 7124 |
| Calculated MW | 17/26 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC=1:50-200

Note: This product is for in vitro research use only