

## CA125 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90050Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CA125 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	250-350/2008
Conjugates	Bio
Target Name	CA125
Other Names	CA125
Accession No.	Uniprot:Q8WXI7Gene ID:94025
Uniprot	Q8WXI7
GeneID	94025
Calculated MW	200
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only