

TNF beta Polyclonal Antibody Biotin Conjugated

Catalog No: #C90051Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF beta Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	91-205/205
Conjugates	Bio
Target Name	TNF beta
Other Names	TNF beta
Accession No.	Uniprot:P01374Gene ID:4049
Uniprot	P01374
GeneID	4049
Calculated MW	18
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only