

alpha 1 Antitrypsin Polyclonal Antibody Biotin Conjugated



Catalog No: #C90053Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	alpha 1 Antitrypsin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/418
Conjugates	Bio
Target Name	alpha 1 Antitrypsin
Other Names	alpha 1 Antitrypsin
Accession No.	Uniprot:P01009Gene ID:5265
Uniprot	P01009
GeneID	5265
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only