

TGF beta 2 Propeptide Polyclonal Antibody Biotin Conjugated



Catalog No: #C90055Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGF beta 2 Propeptide Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	154-197/414
Conjugates	Bio
Target Name	TGF beta 2 Propeptide
Other Names	TGF beta 2 Propeptide
Accession No.	Uniprot:P61812Gene ID:7042
Uniprot	P61812
GeneID	7042
Calculated MW	12/45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only