

# beta Amyloid 1-42 Polyclonal Antibody Biotin Conjugated



Catalog No: #C90058Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	beta Amyloid 1-42 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-42/42
Conjugates	Bio
Target Name	beta Amyloid 1-42
Other Names	beta Amyloid 1-42
Accession No.	Uniprot:P05067Gene ID:351
Uniprot	P05067
GeneID	351
Calculated MW	4.4
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only