

progesterone receptor Polyclonal Antibody Biotin Conjugated

Catalog No: #C90061Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	progesterone receptor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/933
Conjugates	Bio
Target Name	progesterone receptor
Other Names	progesterone receptor
Accession No.	Uniprot:P06401Gene ID:5241
Uniprot	P06401
GeneID	5241
Calculated MW	103
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only