

## ERBB2/HER2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90072Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ERBB2/HER2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu
Immunogen Type	1201-1255/1255
Conjugates	Bio
Target Name	ERBB2/HER2
Other Names	ERBB2/HER2
Accession No.	Uniprot:P04626Gene ID:2064
Uniprot	P04626
GeneID	2064
Calculated MW	185
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only