

SARS S-protein Polyclonal Antibody Biotin Conjugated

Catalog No: #C90075Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SARS S-protein Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	SARS coronavirus
Immunogen Type	401-500/1255
Conjugates	Bio
Target Name	SARS S-protein
Other Names	SARS S-protein
Accession No.	Uniprot:P59594Gene ID:
Uniprot	P59594
Calculated MW	139
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only