

# GTP cyclohydrolase 1 Polyclonal Antibody Biotin Conjugated



Catalog No: #C90077Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GTP cyclohydrolase 1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	34-110/250
Conjugates	Bio
Target Name	GTP cyclohydrolase 1
Other Names	GTP cyclohydrolase 1
Accession No.	Uniprot:P30793Gene ID:2643
Uniprot	P30793
GeneID	2643
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only