

Tubulin-alpha (Loading Control) Polyclonal Antibody FITC Conjugated



Catalog No: #C90083FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tubulin-alpha (Loading Control) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	375-448/448
Conjugates	FITC
Target Name	Tubulin-alpha (Loading Control)
Other Names	Tubulin-alpha (Loading Control)
Accession No.	Uniprot:Q71U36Gene ID:7846
Uniprot	Q71U36
GeneID	7846
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only