

## VEGFR1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90090Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VEGFR1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1162-1260/1333
Conjugates	Bio
Target Name	VEGFR1
Other Names	VEGFR1
Accession No.	Uniprot:P53767Gene ID:14254
Uniprot	P53767
GeneID	14254
Calculated MW	147
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only